

# Disposable Eckardt TwinLight Chandelier

(27 gauge / 0.4 mm)

Article No: 3269.MBD27



**Giuseppe Scarpa,**  
Director of Ophthalmology  
Department, Treviso Hospital,  
Veneto region, Italy.

[Watch surgical video](#)



**“The DORC Eckardt TwinLight Chandelier allows for maximum flexibility in the choice of insertion sites, optimizing the light beam for different vitreoretinal diseases, providing homogeneous illumination and minimizing shadow zones. This product offers an opportunity to regulate light cones by modifying insertion depth of the fibers. Thanks to their flexibility and small diameter, the chandeliers don’t interfere with scleral indentation. It is also suitable for bi-manual management of not only posterior pole diseases, but additionally, most of the other manoeuvres in which it is necessary to have different retinal quadrants constantly under visual control, such as exchanges, aspiration of subretinal fluid and removal of foreign bodies.”**